

TOSHIBA

No. 2D730-124E*D

**SERVICE MANUAL
FOR
DIAGNOSTIC ULTRASOUND SYSTEM**

TOSBEE

**MODEL SSA-240A
(2D730-124E*D)**

TOSHIBA CORPORATION

© TOSHIBA CORPORATION 1990

ALL RIGHTS RESERVED

REVISION RECORD

REV.	DATE (MM/YY)	REASON /AUTHOR	PAGE CHANGED	SER. No.	DOC. PRODUCT.
INI.	07/ '90		-----		KD-M
*A	12/ '90	Mr. Yamazaki			
*B	02/ '92	Mr. Onodera			
*C	06/ '92	Mr. Sato			
*D	08/ '94	Mr. Hirano			KD-WW

To service personnel:

This service manual describes the diagnostic ultrasound system model SSA-240A of the α series. This manual is divided into 2 parts to provide service information for the SSA-240A at the PWB and unit levels.

Part 1: System description volume

Describes the system using block diagrams.

Part 2: Troubleshooting volume

Describes the troubleshooting procedures using the test system.

Precautions concerning PWBs

Like the Y series, SRAM backed up by a battery is used as storage memory on the CPU PWB for initialization data, measurement data, etc. When replacing this PWB, do not place it on a metal surface or the battery may be shorted, resulting in loss of data.

CONTENTS

PART 1 SYSTEM DESCRIPTION VOLUME

	Page
1. OUTLINE -----	1-1-1
2. DISASSEMBLY AND INSTALLATION PROCEDURES -----	1-2-1
2.1 Name of Each Part -----	1-2-1
2.2 Removing the Observation Monitor -----	1-2-2
2.3 Removing the Sheet Keyboard of the Main Panel -----	1-2-3
2.4 Removing the Power Supply Unit (for TAMS) -----	1-2-4
2.5 Removing PWBs (for TAMS) -----	1-2-5
2.6 Layout of PWBs -----	1-2-6
2.7 Removing the ECG AMP PWB (optional) (for TAMS) -----	1-2-7
2.8 Removing the Power Supply Unit (for TMSE) -----	1-2-8
2.9 Removing PWBs (for TMSE) -----	1-2-9
2.10 Removing the ECG AMP PWB (optional) (for TMSE) -----	1-2-10
3. OUTLINE OF EACH PWB -----	1-3-1
3.1 T&R Unit -----	1-3-1
3.1.1 PROBE SELECTOR -----	1-3-2
3.1.2 PULSER -----	1-3-2
3.1.3 MTX SWT -----	1-3-3
3.1.4 R-DELAY -----	1-3-3
3.1.5 RECEIVER -----	1-3-4
3.1.6 MECHA-PULSER -----	1-3-4
3.2 DSC Unit -----	1-3-5
3.2.1 PANEL-IF -----	1-3-5
3.2.2 ECG/LSR-IF -----	1-3-6
3.2.3 DSC-IO/FM PWB (DSC-IO section) -----	1-3-7

CONTENTS - Cont'd

	Page
3.2.4 DSC-IO/FM PWB (DSC-FM section) -----	1-3-9
3.2.5 CPU -----	1-3-10
3.2.6 MECHA-CONT -----	1-3-10
3.2.7 RPG/TR CONT -----	1-3-12
3.2.8 ECG AMP -----	1-3-12
4. BLOCK DIAGRAM -----	1-4-1
4.1 Overall Block Diagram -----	1-4-1
4.1.1 Electronic scan mode -----	1-4-2
4.1.2 Mechanical sector mode -----	1-4-3
4.1.3 LSR M mode -----	1-4-4
4.1.4 VTR mode -----	1-4-5
4.2 Function Block Diagram -----	1-4-6
4.2.1 T&R unit block diagram -----	1-4-6
4.2.2 DSC unit block diagram -----	1-4-7
4.2.3 VIDEO SELECTOR block diagram -----	1-4-8
4.2.4 ECG, LSR block diagram -----	1-4-9
4.2.5 SYSTEM CONTROL block diagram -----	1-4-10
4.3 T & R Unit Block Diagram -----	1-4-11
4.3.1 PROBE SELECTOR block diagram -----	1-4-12
4.3.2 PULSER block diagram -----	1-4-13
4.3.3 MTX SW block diagram -----	1-4-14
4.3.4 R-DELAY block diagram -----	1-4-15
4.3.5 RECEIVER block diagram -----	1-4-16
4.3.6 MECHA-PULSER block diagram -----	1-4-17
4.4 DSC Unit Block Diagram -----	1-4-18
4.4.1 PANEL-IF block diagram -----	1-4-18
4.4.2 ECG/LSR-IF block diagram -----	1-4-20

It's a preview. You can download the full file by clicking the link below.

<https://shopservicemanual.com/>

Service Manuals from 2\$