


1064nm / 532nm Q-Switched Nd:YAG Laser

LUCID^Q-PTPTM

SERVICE TRAINING MANUAL

Contact Information



BISON Medical CO., LTD.

(Gasan-dong, Ace High-end Tower 6), 1801, 1802, 234, Beotkkot-ro, Geumcheon-gu, Seoul, Korea

Tel : 82-2-865-7121, Fax : 82-2-865-7131

bison@bisonmedical.com

www.bisonmedical.com

Your Dealer



Document NO : BS-SM-QP-16

Revision : 1.000

Date : AUG.24.2015



© 2010 BISON Medical CO. LTD.

All rights reserved. Nothing from this publication may be copied, stored in an automated data file, or made public, in any form or in any way, be it electronically, mechanically, by photocopying, recordings or in any other way, without prior written permission from BISON Medical CO., LTD.

Introduction

This content can be seen by the authorized Service Engineer of BISON Medical, and unlawful copy or distribution is not permitted.

This Service Training Manual is composed of total 6 sections. All Service Engineers should be fully trained to be used at the term, definition, and contents. Correct understanding of the contents enables the appropriate preparation when A/S is occurred. This preparation enables to reduce the Service time.

It is recommended to read the User Manual first. In User Manual, cautions in the use of equipment including the installation method and operation method of the equipment etc. are explained in detail.

Please inquire to BISON Medical if you have any questions or insufficient parts on the Training Manual.

Trademark/copyright information

BISON Medical and LUCID Q-PTP registered trademarks of BISON Medical CO.,LTD.

All other trademarks or registered trademarks are the property of their respectively owners



Table of Contents

Table of Contents	IV
1. System Introduction	1
2. Safety Guideline.....	2
2.1. Danger of burn	2
2.2. Protection of a patient.....	2
2.3. Danger of explosion and fire.....	3
2.4. Danger of lethal voltage	3
2.5. Danger of exposure of laser light to the eye.....	3
2.6. Protective Earth of the equipment.....	4
2.7. Confirmation of warnings.....	4
2.8. Disposal	4
2.9. General precautions while in operating.....	5
2.10. Precautions while operation.....	6
2.11. Safety system	7
2.11.1. Key Switch	7
2.11.2. Emergency Switch.....	7
2.11.3. Circuit Breaker	7
2.11.4. Remote Interlock.....	7
2.11.5. Errors.....	8
2.12. Labels	9
3. System Description	10
3.1. Specifications.....	10
3.2. Laser system	10
3.3. Aiming Beam specification	11

- 3.4. Environmental condition..... 11**
- 3.5. Exterior view 12**
- 3.6. System Components 13**
- 3.7. System Diagram..... 14**
- 3.8. System Boards..... 15**
 - 3.8.1. Main Board..... 15
 - 3.8.2. Beam Shutter Board..... 17
 - 3.8.3. Q-YAG Board..... 17
 - 3.8.4. Power Board 18
 - 3.8.5. Jog Shuttle Board..... 19
 - 3.8.6. Relay Board..... 20
 - 3.8.7. PD Board 21
- 3.9. Laser Head 22**
 - 3.9.1. Specifications..... 22
 - 3.9.2. Laser Head Layout..... 22
 - 3.9.3. Wavelength Control 23
 - 3.9.4. Head Connections 23
 - 3.9.5. Plug Description 24
- 3.10. Arm Socket..... 26**
- 3.11. Touch LCD 27**
- 3.12. Water Pump..... 27**
- 3.13. High voltage power supply(HVPS) 28**
- 3.14. Chiller 29**
- 3.15. Flow Switch..... 30**
- 3.16. Temperature Controller..... 31**
- 3.17. Rear panel connections..... 32**

It's a preview. You can download the full file by clicking the link below.

<https://shopservicemanual.com/>

Service Manuals from 2\$